

# EMT Nozzle

## Enhanced Mist Technology



Highly innovative and patented 2 piece insert design. Specifically designed for low VOC and water-based formulations using compressed air systems.

Improved performance and control. For personal and home care applications.

### Currently available with:

- ✔ Lindal standard actuators:
  - Neo** - for male and female valves
  - ST370** - for male valves
  - 370** - for female valves
- ✔ Lindal standard valves:
  - LI** - male valve
  - RT** - female valve
  - BOV** - Bag on Valve (male & female)
  - LIB** - Ball Valve (male valve)



### More info:

Visit: [www.lindalgroup.com](http://www.lindalgroup.com)

Email: [products@lindalgroup.com](mailto:products@lindalgroup.com)

Contact: Your local sales consultant



Your Innovation Partner

# EMT Nozzle

Technical details  
Enhanced Mist Technology

## Available configurations

EMT300





EMT400

## Technical details

T130 – Lindal Standard non-Mist nozzle

EMT300

EMT400



**i** All performance tests done with demineralized water – Spray diameter test done with a spraying distance of 15cm. Tests results are informative only and are non-binding.

### More info:

Visit: [www.lindalgroup.com](http://www.lindalgroup.com)

Email: [products@lindalgroup.com](mailto:products@lindalgroup.com)

Contact: Your local sales consultant



Your Innovation Partner